DX2 REM550 (DYNAMIC CONTROLS)

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
Afbeelding 8.2

English

© 2017 Scout Mobility B.V. All rights reserved.

The information provided herein may not be reproduced and/or published in any form, by print, photoprint, microfilm or any other means whatsoever (electronically or mechanically) without the prior written authorization of Scout Mobility B.V..

The information provided is based on general data concerning the constructions known at the time of the publication of this manual. Scout Mobility B.V. executes a policy of continuous improvement and reserves the right to changes and modifications.

The information provided is valid for the product in its standard version. Scout Mobility B.V. cannot be held liable for possible damage resulting from specifications of the product deviating from the standard configuration.

The available information has been prepared with all possible diligence, but Scout Mobility B.V. cannot be held liable for possible errors in the information or the consequences thereof. Scout Mobility B.V. accepts no liability for loss resulting from work executed by third parties.

Names, trade names, etc. used by Scout Mobility B.V. may not, as per the legislation concerning the protection of trade names, be considered as being available.

2017-01

Table of content

Preface	17
1. The controller	18
1.1 R-net controller	18
1.2 Display screen	20
2. Driving the wheelchair with the controller	21
2.1 Switching the controller on or off	21
2.2 Driving the wheelchair	21
2.3 Speed.....	21
2.4 Driving Profile.....	21
2.5 Selecting a profile.....	21
3. Operating the powered adjustment options	22
3.1 Selecting the desired adjustment option	22
4. Lights	23
5. Troubleshooting	23
6. Locking the controller	23
7. Technical specifications	24
8. Technical diagrams	25
8.1 Technical Diagram	25
8.2 Technical diagram battery charger.....	25

Preface

This manual

This user manual describes the use of the power wheelchair via the controller. The total user manual for this power wheelchair consists of three booklets. Read the entire user manual carefully before taking the product into use. The information provided by this manual is essential for the safe use and proper care (cleaning) of the wheelchair.

- The general usermanual
- The user manual for the seating system (AGILO)
- The user manual for the controller (this booklet)

If one of these user manuals was not included with your wheelchair, please contact your dealer immediately. In addition to this user manual, there is also a service manual for qualified specialists.

This general user manual will refer, where necessary, to one of the other user manuals. This is indicated as follows:

GENERAL: Refers to the general user manual.

CONTROLLER: Refers to the user manual for the controllers.

AGILO: Refers to the user manual for the seating system.

1. The controller

1.1 R-net controller

A controller will usually have three basic functions:

- Driving and steering a wheelchair
- Operating powered seat adjustments
- Charging the wheelchairs batteries

There are many different control systems for wheelchairs on the market. If the controller on your wheelchair does not resemble the one in figure 1.1, contact your dealer.


Figure 1.1

DX2 of Dynamic Controls is a collective name for the entire control system of the wheelchair. The wheelchair is operated by means of a controller, which includes the following components:


Figure 1.2

Part	Function
A. Joystick	In the 'drive' mode: driving and steering In the 'adjustment options' mode:: - Left/right to select the adjustment options - Front/back to select the adjustment options
B. On/off switch	Switching the controller on or off
C. Horn	Warning signal with sound
D. LCD color display screen	Display and feedback
E. "Adjustments" button	To select powered options. The display shows the chosen powered adjustment option.
F. "Function" button	Depending the mode you are working in - In 'drive' mode; slower / faster - In 'adjustment options' mode; selection of the powered adjustment option.
G. Charge connector	Input for the battery charger
H. Charge plug for the battery charger	Connector for the battery charger
I. Lights button	Push >2sec: Switching the lights on or off
J. Hazard lights	Push >2sec: Switching on or off hazard lights
K1: Direction indicator left	Switches the left direction indicator on or off
K2: Direction indicator right	Switches the right direction indicator on or off
L. "Profile" button	Select driving profile

1.2 Display screen


Figure 1.3: Display

Part	Function
D1: Battery indicator	Displays the power level of the battery
D2: Active direction indicators	Displays if the direction indicators/hazard lights are on
D3: Active side lights	Displays if the side lights are on
D4: Fault indicator	Displays a code number if a fault occurs
D5: Clock	Displays the time
D6: Speed indicator	Displays a graph of the actual speed
D7: Profile	Displays the selected profile
D8: Maximum speed	Displays the maximum speed limit as set by the user
E: Adjustment options	Displays the chosen powered adjustment option
F: Functions	Displays the selectable functions. Depending the mode you are working in.
L1: Next profile	To select the next drive profile
L2: previous profile	To select the previous drive profile

2. Driving the wheelchair with the controller

2.1 Switching the controller on or off

To be able to drive or operate the powered adjustment options of the wheelchair, the controller must be switched on. Press the on/off button (B in figure 1.2).

2.2 Driving the wheelchair

Driving an power wheelchair is done by operating a joystick. Move the joystick forwards and the wheelchair will also move forwards. Steer left and right and the wheelchair will turn.

2.3 Speed

The maximum speed can be controlled by the speed regulator on the controller (L in figure 1.2). The speed bar on the display screen will display the maximum speed (D8 in figure 1.3). Speed can be controlled with the joystick while driving. If the joystick is moved a little, the wheelchair will move more slowly.

NB: The maximum speed can vary by each drive profile!

2.4 Driving Profile

This controller can also be set to make the wheelchair suitable for different driving profiles or environments. For example, selecting the profile for indoors will ensure that the wheelchair will react more 'calmly'. Once outside, the profile can be adjusted to a more 'robust' setting. The name (number) of the current profile can be found in the display (D7 in figure 1.3). DX2 REM550 allows you to adjust the maximum speed within different profiles or environments.

2.5 Selecting a profile

To switch to a different profile, just press the 'profile' button (L in figure 1.2). Press the "plus" or "minus" till you have reached the desired drive profile.

3. Operating the powered adjustment options

Not every wheelchair has been equipped with powered adjustment options. We make a distinction between four different adjustment options to the seating system:

- Leg rest adjustment
- Powered seat lift adjustment
- Tilt adjustment
- Backrest adjustment


Figure 3.1 powered adjustment options

3.1 Selecting the desired adjustment option

1. Switch on the controller
2. Press on the 'adjustment' button to select the 'adjustment options' mode (E in figure 1.2).
3. Move the joystick to the left or the right to select the desired adjustment option (or use the 'function button' (F in figure 1.2) . The selected adjustment option will be visible on the display screen.
4. Moving the joystick forwards and/or backwards will activate the selected powered adjustment option (see table 3.1). Move the joystick forwards or backwards until the desired adjustment option has been attained.

	Move joystick backwards	Move joystick forwards
Tilt adjustment	The entire chair will tilt backwards	The entire chair will tilt forwards
Seat lift adjustments	The entire chair will be raised	The entire chair will be lowered
Backrest adjustments	The backrest will tilt backwards	The backrest will tilt forwards
Legrest adjustments	The legrest angle will increase, the footplate will be raised	The legrest angle will decrease, the footplate will be lowered

Table 3.1 powered adjustments

5. To return to the 'drive' mode: Press on the 'Profile' button to select the 'drive' mode.

Note: If you move the chair while adjusting the powered seatlift option, the speed will be reduced.

4. Lights

Not every wheelchair has been equipped with lights.

- Lights (I in figure 1.1)
- Hazard Lights (J in figure 1.1)
- Direction indicators (K1 & K2 in figure 1.1)

5. Troubleshooting

If the wheelchair will not function while the batteries are fully charged, check the following points before consulting your dealer:

- Switch the controller off and then switch it on again. Check to see if the malfunction has been solved.
- Check if the free wheel switch was switched to Drive.
- Check if the joystick was in the 0 position when the controller was switched on. In other words, the joystick must not be moved when the controller is being switched on or off.

Malfunctions list

An extensive list of possible malfunctions can be found in the appendix of the service manual for qualified specialists (also available on www.scoutmobility.nl)

6. Locking the controller

To lock the controller DX2 REM550

- Put the controller on (B in fig 6.1).
- Push >4 sec the on/off button (B in fig 6.1)

To unlock the controller DX2 REM550

- Put the controller on (B in fig 6.1)
- Push 2x the claxon (C in fig 6.1)


Figure 6.1

7. Technical specifications

Technical specifications DX2 REM 550				
Parameter	Value			
Operating Voltage range	18V -32V DC (nom. 24V)			
Charger rating	Max 12 A RMS Continuous, limited by DXBUS Rating			
Quiescent Current	< 0.25mA Off, typically 200mA On			
Protection Rating	IPx4			
Parameter	Min	Nominal	Max	Units
Tube mount diameter	15 (5/8)	19 (3/4)	22 (7/8)	mm (in)
Operating Temperature Range	-25 (-13)		50 (122)	°C (°F)
Storage Temperature Range	-40 (-40)		70 (158)	°C (°F)
Operating Humidity Range	0		95	%RH

Technical specifications DX2 Power Module				
Parameter	Min	Nominal	Max	Units
Battery Voltage	18	24	32	V
Battery Saver				
- High threshold		21		
- Low threshold		18		
Quiescent Power (relay de-energised)		2	5	W
Standby Power		15	20	mW
Thermal limit		70		°C
DX BUS				
CANH and CANL tolerated Voltage	0		38	V
DX BUS supply current @ 40°C	12	14.9		A
Motor outputs				
Motor Type		24V		DC
Maximum output voltage	23.5			V
Continuous motor current (@ 20°C ambient temperature)				
- DX2-PMA70/DX2-PMA70L		20		A
- DX2-PMA90/DX2-PMA90L		30		A
Peak motor current				
- DX2-PMA70/DX2-PMA70L	30*	70**	75***	A
- DX2-PMA90/DX2-PMA90L	35*	90**	120***	A
* 15 min @ 20°C initial case temperature				
** 30 sec @ 20°C initial case temperature				
*** 0.1-1.5 sec (programmable boost time)				
PWM frequency	16	19.6	20	kHz
Stall timer (programmable)	2	15	60	s
Parkbrake outputs				
Output voltage	$V_{bat} \cdot 1.1$		V_{bat}	V
Output current	0.6	1		A
Configuration: single/dual				
Lighting outputs (70L and 90L variants only)				
Supply Voltage	18	24	32	V
Output Power (per output)			50	W
Output Current (per output)			2	A
On Board charger supply				
Voltage		24		V
Current			15	A
Inhibit	Via DCI input			

8. Technical diagrams

8.1 Technical Diagram

The technical diagram can also be found in the battery box.

8.2 Technical diagram battery charger

The controller's standard configuration includes a '3-pin connection'. Ensure that the battery charger is properly connected so that the 'negative pole' and the 'inhibit' are connected, enabling the system to prevent the wheelchair from moving when the battery is being charged.


Figure 8.2

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If there is no information here you can contact Scout Mobility B.V. for the nearest dealer address.