
SHARK (DYNAMIC CONTROLS)

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
Afbeelding 8.2

English

© 2017 Scout Mobility B.V. All rights reserved.

The information provided herein may not be reproduced and/or published in any form, by print, photoprint, microfilm or any other means whatsoever (electronically or mechanically) without the prior written authorisation of Scout Mobility B.V..

The information provided is based on general data concerning the constructions known at the time of the publication of this manual. Scout Mobility B.V. executes a policy of continuous improvement and reserves the right to changes and modifications.

The information provided is valid for the product in its standard version. Scout Mobility B.V. cannot be held liable for possible damage resulting from specifications of the product deviating from the standard configuration.

The available information has been prepared with all possible diligence, but Scout Mobility B.V. cannot be held liable for possible errors in the information or the consequences thereof. Scout Mobility B.V. accepts no liability for loss resulting from work executed by third parties.

Names, trade names, etc. used by Scout Mobility B.V. may not, as per the legislation concerning the protection of trade names, be considered as being available.

2017-01

Preface.....	16
1. The controller	17
1.1 Shark.....	17
2. Driving the wheelchair with the controller	19
2.1 Switching the controller on or off	19
2.2 Driving the wheelchair.....	19
2.3 Speed.....	19
3. Operating the powered adjustment options	20
3.1 Selecting the desired adjustment option	20
4. Lights	21
5. Troubleshooting.....	21
6. Locking the controller	21
7. Technical specifications.....	21
8. Technical diagrams	23
8.1 Technical Diagram	23
8.2 Technical diagram battery charger.....	23

Preface

This manual

This user manual describes the use of the power wheelchair via the controller. The total user manual for this power wheelchair consists of three booklets. Read the entire user manual carefully before taking the product into use. The information provided by this manual is essential for the safe use and proper care (cleaning) of the wheelchair.

- The general usermanual
- The user manual for the seating system (AGILO)
- The user manual for the controller (this booklet)

If one of these user manuals was not included with your wheelchair, please contact your dealer immediately. In addition to this user manual, there is also a service manual for qualified specialists.

This general user manual will refer, where necessary, to one of the other user manuals. This is indicated as follows:

GENERAL: Refers to the general user manual.

CONTROLLER: Refers to the user manual for the controllers.

AGILO: Refers to the user manual for the seating system.

1. The controller

1.1 Shark

A controller will usually have three basic functions:

- Driving and steering a wheelchair
- Operating powered seat adjustments
- Charging the wheelchairs batteries

There are many different control systems for wheelchairs on the market. If the controller on your wheelchair does not resemble the one in figure 1.1, contact your dealer.


Figure 1.1

Shark of Dynamic Controls is a collective name for the entire control system of the wheelchair. The wheelchair is operated by means of a controller, which includes the following components:


Figure 1.2 Left REM421; Right REM420

Part	Function
A: Joystick	In the 'drive' mode: driving and steering In the 'adjustment options' mode:: - Left/right to select the adjustment options - Front/back to select the adjustment options
B: On/off switch	Switching the controller on or off
C: Horn	Warning signal with sound
D1: Battery indicator	Displays the power level of the battery
D2: Maximum speed	Displays the maximum speed limit as set by user.
D3: Indicator / hazard lights display	Displays status of the indicators and hazard lights (on/off)
D4: Fault indicator	Displays fault code in case of system error
E: "Adjustments" button	To select the powered seating adjustment (1, 2, 3, 4).
F1: Speed regulator	Reduce driving speed (slower)
F2: Speed regulator	Increase driving speed (faster)
G: Charge connector	Input for the battery charger
H: Charge plug for the battery charger	Connector for the battery charger
I: Lights button	Switching the lights on/off
J: Hazard lights	Warning signal with lights
K1: Direction indicator left	Switches the left direction indicator on/off
K2: Direction indicator right	Switches the right direction indicator on/off

2. Driving the wheelchair with the controller

2.1 Switching the controller on or off

To be able to drive or operate the powered adjustment options of the wheelchair, the controller must be switched on. Press the on/off button (B in figure 1.2).

2.2 Driving the wheelchair

Driving a power wheelchair is done by operating a joystick. Move the joystick forwards and the wheelchair will also move forwards. Steer left and right and the wheelchair will turn.

2.3 Speed

The maximum speed can be controlled by the 'speed regulator' on the controller (F in figure 1.2). Speed can be controlled with the joystick while driving. If the joystick is moved a little, the wheelchair will move more slowly.


Figure 2.1 Driving profiles

3. Operating the powered adjustment options

Not every wheelchair has been equipped with powered adjustment options. We make a distinction between four different adjustment options to the seating system:

- Legrest adjustment
- Powered seatlift adjustment
- Tilt adjustment
- Backrest adjustment

Figure 3.1 Powered adjustment options

3.1 Selecting the desired adjustment option

1. Switch on the controller
2. Press on the ‘adjustment’ button to select the ‘adjustment options’ mode (E in figure 1.2).
3. Push the “1” or the “2” to select the desired adjustment. The corresponding number enlightens on the display.
4. Moving the joystick forwards and/or backwards will activate the selected powered adjustment option (see table 3.1). Move the joystick forwards or backwards until the desired adjustment option has been attained.

	Move joystick backwards	Move joystick forwards
Tilt adjustment	The entire chair will tilt backwards	The entire chair will tilt forwards
Seat lift adjustments	The entire chair will be raised	The entire chair will be lowered
Backrest adjustments	The backrest will tilt backwards	The backrest will tilt forwards
Legrest adjustments	The legrest angle will increase, the footplate will be raised	The legrest angle will decrease, the footplate will be lowered

Table 3.1 Powered adjustments

5. To return to the ‘drive’ mode: Press on the ‘Speed Regulator’ button to select the ‘drive’ mode.

Note: If you move the chair while adjusting the powered seatlift option, the speed will be reduced.

4. Lights

Not every wheelchair has been equipped with lights.

- Lights (I in figure 1.1)
- Hazard Lights (J in figure 1.1)
- Direction indicators (K1 & K2 in figure 1.1)

5. Troubleshooting

If the wheelchair will not function while the batteries are fully charged, check the following points before consulting your dealer:

- Switch the controller off and then switch it on again. Check to see if the malfunction has been solved.
- Check if the free wheel switch was switched to Drive.
- Check if the joystick was in the 0 position when the controller was switched on. In other words, the joystick must not be moved when the controller is being switched on or off.

Malfunctions list

An extensive list of possible malfunctions can be found in the appendix of the service manual for qualified specialists (also available on www.scoutmobility.nl)

6. Locking the controller

To lock the controller Shark

- Push >4 sec the on/off button (B in fig 6.1) while the chair is turned on.
- The DX2 system turns off immediately when the on/off button is pushed.
- After 4 seconds the 'maximum speed' (D2 in fig 1.2) and the Lock LED turn on for 1 second to indicate that the DX2 system is now locked.

To unlock the controller DX2 REM420/421

- Press the on/off button
- Press the horn button twice (within 10 sec).
- The wheelchair is now unlocked again and can be used again.

7. Technical specifications

Technical specifications Shark				
Parameter	Value			
Operating Voltage range	18V -32V DC (nom. 24V)			
Charger rating	Max 12 A RMS Continuous, limited by Shark BUS Rating			
Quiescent Current	< 0.25mA Off, typically 200mA On			
Protection Rating	IPx4			
Operating force	Joystick < 2.2N, Buttons < 3N			
Parameter	Min	Nominal	Max	Units
Tube mount diameter	15 (5/8)	19 (3/4)	22 (7/8)	mm (in)
Operating Temperature Range	-25 (-13)		50 (122)	°C (°F)
Storage Temperature Range	-40 (-40)		65 (149)	°C (°F)
Operating Humidity Range	0		90	%RH

Technical specifications SHARK Power Module				
Parameter	Min	Nominal	Max	Units
Compatible Battery Supply	24V supply, 2 x 12V in series, circuit breaker protected, gel cell type of recommended minimum capacity 20 Ah.			
Compatible Motor	24V DC permanent magnet type, typically rated 100-300 W.			
Compatible Park Brake	Either 2 x 24V connected in parallel or 2 x 12V connected in series.			
Shark operating voltage (V_{bat})	18	24	32	V
Reverse Supply Voltage	-32			V
Continuous current (@ 20°C ambient)				
DK-PMB	16	18		A
Peak current (<60 sec @ 20°C initial)				
DK-PMB	55	60	63	A
Boost current (< 10 sec on top of peak current)				
DK-PMB	0		15	A
Parkbrake outputs				
Output voltage	$V_{\text{bat}} 1.2$		V_{bat}	V
Output current	0.6	0.7		A
Quiescent current (idle)		120	150	mA
Drive Control Input				
Short Circuit Output Current	8		15	mA
Maximum Input Voltage Range	0		50	V
SHARK Communications Bus				
SBHi, SBLo Voltage Range	0		38	V
Max Output Current (SBB+)		8	12 (RMS)	A
Controller Battery Charger Max Current		8	12 (RMS)	A
Power Module OBC Max Current			5 (RMS)	A
Seat 1 / 2 Current (60 second rating, programmable)	0	6	12	A
Seat 1 / 2 Output Voltage (with Seat X Speed set to 100 %, proportional)		Battery Voltage - 1.0V		V
Nominal Left/Right/Head Light Bulb Voltage		24		V
Lighting Circuit Supply Voltage (Shark Port Pins 3 + 12)		Battery Voltage - 1.0V		V
Lighting Circuit Supply Current (Shark Port Pins 3 + 12)			8	A
Headlight Output Current Rating		3		A
Left/Right Indicator Output Current Rating		3		A

8. Technical diagrams

8.1 Technical Diagram

The technical diagram can also be found in the battery box.

8.2 Technical diagram battery charger

The controller's standard configuration includes a '3-pin connection'. Ensure that the battery charger is properly connected so that the 'negative pole' and the 'inhibit' are connected, enabling the system to prevent the wheelchair from moving when the battery is being charged.


Figure 8.2

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Technical specifications Shark				
Parameter	Value			
Operating Voltage range	18V -32V DC (nom. 24V)			
Charger rating	Max 12 A RMS Continuous, limited by Shark BUS Rating			
Quiescent Current	< 0.25mA Off, typically 200mA On			
Protection Rating	IPx4			
Operating force	Joystick < 2.2N, Buttons < 3N			
Parameter	Min	Nominal	Max	Units

Tube mount diameter	15 (5/8)	19 (3/4)	22 (7/8)	mm (in)
Operating Temperature Range	-25 (-13)		50 (122)	°C (°F)
Storage Temperature Range	-40 (-40)		65 (149)	°C (°F)
Operating Humidity Range	0		90	%RH

Technical specifications SHARK Power Module				
Parameter	Min	Nominal	Max	Units
Compatible Battery Supply	24V supply, 2 x 12V in series, circuit breaker protected, gel cell type of recommended minimum capacity 20 Ah.			
Compatible Motor	24V DC permanent magnet type, typically rated 100-300 W.			
Compatible Park Brake	Either 2 x 24V connected in parallel or 2 x 12V connected in series.			
Shark operating voltage (V_{bat})	18	24	32	V
Reverse Supply Voltage	-32			V
Continuous current (@ 20°C ambient)				
DK-PMB	16	18		A
Peak current (<60 sec @ 20°C initial)				
DK-PMB	55	60	63	A
Boost current (< 10 sec on top of peak current)				
DK-PMB	0		15	A
Parkbrake outputs				
Output voltage	$V_{bat} 1.2$		V_{bat}	V
Output current	0.6	0.7		A
Quiescent current (idle)		120	150	mA
Drive Control Input				
Short Circuit Output Current	8		15	mA
Maximum Input Voltage Range	0		50	V
SHARK Communications Bus				
SBHi, SBLo Voltage Range	0		38	V
Max Output Current (SBB+)		8	12 (RMS)	A
Controller Battery Charger Max Current		8	12 (RMS)	A
Power Module OBC Max Current			5 (RMS)	A
Seat 1 / 2 Current (60 second rating, programmable)	0	6	12	A
Seat 1 / 2 Output Voltage (with Seat X Speed set to 100 %, proportional)		Battery Voltage - 1.0V		V
Nominal Left/Right/Head Light Bulb Voltage		24		V
Lighting Circuit Supply Voltage (Shark Port Pins 3 + 12)		Battery Voltage - 1.0V		V
Lighting Circuit Supply Current (Shark Port Pins 3 + 12)			8	A
Headlight Output Current Rating		3		A
Left/Right Indicator Output Current Rating		3		A

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
Abbildung 8.2

Dealer contact details:

If there is no information here you can contact Scout Mobility B.V. for the nearest dealer address.

Scout Mobility B.V.
Kanaalstraat 101D
5711 EG Someren
The Netherlands

T +31(0) 6 44500578
info@scoutmobility.nl
www.scoutmobility.nl